



preOx[®]-HUT

Rapid Liquid Urease Test

HELICOBACTER PYLORI

preOx[®]-HUT for detection of Helicobacter pylori infection in stomach is based on the activity of the enzyme urease in a biopsy specimen.

The biopsy specimen taken from the stomach is tested immediately. The change of color in the test liquid from **yellow** to **red/pink** indicates within 2 to 30 minutes the presence of Helicobacter pylori.



preOx[®]-HUT Advantages

- easy and safe handling
- quick response 2-5 minutes
- cost effective

Order Information

REF. No. 200.100

REF. No. 200.050

Description: H. pylori Liquid Quick Test